

iPhone Kitting Instructions

Step-by-step setup, labeling, SIM, quality check and packing — for CDW Configuration Services

VERSION

1.0 — DRAFT

LAST UPDATED

September 11, 2026

OWNER

BCH CMDI (Clinical Mobile Devices)

DEVICE OS

TBD: iOS version

DRAFT — not yet approved for use. Items highlighted in yellow are open questions BCH must answer before this guide is sent to CDW. Sensitive values (Wi-Fi password, device passcode, ship-to address) are **never** written on this page — BCH provides them to CDW separately.

How this works

Every iPhone on this order is already registered to Boston Children's Hospital in **Apple Business Manager (ABM)**. The first time a phone is powered on and connected to Wi-Fi, it automatically downloads BCH's management profile (MDM) and then installs BCH's apps on its own. Your job is to label each phone, give it cellular service, walk it through setup, set a few settings, open two apps once, check it, and pack it.

0 Prepare

Supplies, Wi-Fi, charge

1 Label

Serial label on the phone's back

2 SIM / eSIM

Cellular service + record numbers

3 Setup Assistant

Power on, enroll in BCH management

4 Hubb app

Open and allow permissions

5 Apps installed

Wait until every app finishes

6 Settings

AutoFill, Live Activities, 5G, App Store

7 Rover & Haiku

Open once, allow prompts

8 Quality check

Verify and log each phone

9 Case & pack

Case on, box, ship

Estimated hands-on time per phone: TBD: minutes . Menu names below match iOS TBD ; they can differ slightly on other versions.

The rules that matter most

If you remember nothing else, remember these. Breaking one usually means the phone has to be erased and set up again.

- 1 Never sign in with an Apple Account (Apple ID)** — not during setup, not in the App Store, not in Settings.
- 2 Never skip Remote Management.** If the *Remote Management* screen does not appear during setup, stop working on that phone and follow Troubleshooting.
- 3 Never sign in to any app** (Hubb, Rover, Haiku, Locker, or any other) with any username or password. Only accept or allow the prompts described here.
- 4 Keep phones on Wi-Fi and on power** until every app has finished installing.
- 5 Do not install, delete or move apps, and do not change settings** that aren't in this guide.
- 6 When something doesn't match this guide, set the phone aside and ask BCH** — don't improvise.

0 Before you start

Supplies

- iPhones — model **TBD: model / storage / color**, quantity **TBD**
- Labels + printer — **TBD: label size / stock** (see Step 1)
- SIM cards or eSIM activation info (see Step 2) and a SIM eject tool if physical SIMs are used
- Cases — **TBD: case make/model**
- Screen protectors — **TBD: yes/no, model**
- Charging cables — **TBD: type (USB-C to USB-C?) and qty**; power adapters **TBD: included?**
- Isopropyl alcohol wipes (to clean the phone back before labeling)
- The device manifest spreadsheet (see Manifest)

From BCH (sent to you separately)

- Wi-Fi network name and password
- Device passcode, if one is required (**TBD: is one required?**)
- SIM assignment sheet / eSIM activation details
- Ship-to address and receiving contact
- Blank manifest template

Wi-Fi requirements

- Internet access with **no sign-in / splash page** (captive portal) — enrollment fails behind one.
- Apple services not blocked: TCP 80, 443, 5223 and 2197 to Apple's network (17.0.0.0/8).

- Enough bandwidth for many phones downloading apps at once (about **TBD** GB per phone).

Get ready

1. ☐ Confirm the serial numbers on the boxes match the order list from BCH. Report any serial that isn't on the list **before** powering it on — it may not be registered to BCH.
2. ☐ Charge phones to at least **50%** before starting, and keep them on power during setup.
3. ☐ Open the manifest and fill in one row per phone as you go.

1 Label the phone

The serial number label goes on the **back of the iPhone itself — not on the case.**

1. ☐ Print the label. It must show:
 - The iPhone serial number as a **Code 128 barcode** **TBD: confirm barcode type** and as readable text
 - **TBD: anything else? (asset tag, "Property of Boston Children's Hospital", phone number, IMEI)**

Label size: **TBD**. The serial number is printed on the iPhone box.

2. ☐ Wipe the back of the phone with an alcohol wipe and let it dry.
3. ☐ Apply the label: **TBD: exact spot — suggested: lower third of the back, centered, horizontal**.
 - Do **not** cover the cameras, flash, or microphone.
 - Keep it flat, straight, with no bubbles or lifted corners. Press firmly for a few seconds.
4. ☐ Make sure you labeled the right phone: after setup, the serial in **Settings** > **General** > **About** must match the label (checked again in Quality check).

2 SIM / eSIM

READ FIRST

iPhone 14 and newer models sold in the US have **no physical SIM tray — they are eSIM only**.
Use Option A or B below based on what BCH confirms: **TBD: physical SIM or eSIM? Carrier?**

Option A — Physical SIM card

1. ☐ With the phone **powered off**, eject the SIM tray with the eject tool.
2. ☐ Insert the SIM card assigned to this phone: **TBD: any SIM in any phone, or a specific SIM per serial?**
3. ☐ Close the tray fully.
4. ☐ Record the SIM's **ICCID** (printed on the SIM card carrier) against the phone's serial in the manifest.

Option B — eSIM

1. ☐ Activation method: **TBD: carrier pre-assigns eSIM to each IMEI/EID (activates on its own during setup) OR a QR code per phone**
2. ☐ During setup, if a *Cellular Setup* or *Carrier eSIM Available* screen appears, tap **Continue** / **Activate**. If using QR codes, scan the code assigned to this phone.
3. ☐ Record the phone's **EID** and **IMEI** (**Settings** > **General** > **About**) in the manifest.

In both cases, record the **phone number** once service is active (**Settings** > **Cellular** → *My Number*) — checked in Quality check.

3 Power on & Setup Assistant

Press and hold the side button until the Apple logo appears. Go through the screens below. **Some screens may be skipped automatically** by BCH's configuration — that's normal. The order can vary slightly.

SCREEN	WHAT TO DO
Hello	Swipe up from the bottom of the screen.
Language	Select English .
Country or Region	Select United States .
Quick Start	Tap Set Up Without Another Device .
Languages / Accessibility / Appearance	Leave the defaults and tap Continue .
Choose a Wi-Fi Network	Select the Wi-Fi network BCH gave you and enter the password.
Activating iPhone / Cellular Setup	Wait. If an eSIM screen appears, see Step 2, Option B.
Data & Privacy	Tap Continue .
Remote Management	Required. Confirm the screen names Boston Children's Hospital TBD: confirm the exact organization name shown , then tap Enroll This iPhone (may say Continue). Wait for it to finish. If this screen never appears, stop — see Troubleshooting.
Create a Passcode	TBD: use the passcode BCH provided, OR tap Passcode Options → Don't Use Passcode
Face ID	Tap Set Up Later .

SCREEN	WHAT TO DO
Transfer Your Apps & Data	Tap Don't Transfer Anything / Set Up as New .
Apple Account	Do not sign in. Tap <i>Forgot password or don't have an Apple Account?</i> → Set Up Later in Settings → Don't Use .
Location Services	Tap Enable Location Services .
Update iPhone Automatically	TBD: Continue (automatic updates on) OR Only Download / Don't install
Siri	Tap Set Up Later .
Screen Time	Tap Set Up Later .
iPhone Analytics / App Analytics	Tap Don't Share .
Display, True Tone, Action Button, Camera Control, Emergency SOS	Leave the defaults and tap Continue .
Welcome to iPhone	Swipe up to reach the Home Screen.
Any screen not listed	Accept the default and tap Continue — unless it asks you to sign in, share data, or skip Remote Management. For those, choose <i>Not Now</i> / <i>Set Up Later</i> / <i>Don't Share</i> .

✓ **Expected:** The Home Screen appears and BCH's apps start installing on their own (gray icons that say *Waiting...* or *Loading...*).

4 Hubb app

Wait until the Hubb icon finishes installing (no gray icon or progress ring), then:

1. ☐ Tap **Hubb** to open it. TBD: does Hubb ask for a login, site code or QR code? If so, what should CDW do?

2. ☐ **Data & Privacy:** tap **Accept**.

3. ☐ **Location:** iOS first offers only *Allow Once / Allow While Using App / Don't Allow*. Tap **Allow While Using App**. When a second prompt offers **Change to Always Allow**, tap it.

4. ☐ For **every other pop-up** (Notifications, Bluetooth, Motion & Fitness, Local Network, etc.) tap **Allow** / **Yes** / **Always Allow**.

5. ☐ Go through the remaining in-app screens and tap **I Understand** on each.

6. ☐ Double-check location: **Settings** > **Apps** > **Hubb** > **Location** must be **Always**, and **Precise Location** must be **on**. Fix it here if the in-app prompt was missed.
✓ **Expected:** Hubb location = Always, Precise Location on.

5 Confirm all apps are installed

Keep the phone on Wi-Fi and power. Apps install on their own — you do not need the App Store.

1. ☐ Check that every app on this list is on the phone: TBD: full list — e.g. Hubb, Rover, Haiku, Locker, CMDI GCT, ...

2. ☐ No icon is gray, shows a progress ring, or says *Waiting...* / *Loading...*

3. ☐ If apps are still waiting after **30 minutes** TBD: confirm time, see Troubleshooting.

6 Device settings

SHORTCUT: THE CMDI GCT APP

The **CMDI GCT** app (installed with the other BCH apps) walks you through these settings screen by screen. Open it and follow its prompts, then use the list below to confirm nothing was missed. **TBD: confirm GCT installs automatically; which steps it covers; add a screenshot**

6.1 AutoFill — use Locker

1. ☐ Go to **Settings** > **General** > **AutoFill & Passwords** .
2. ☐ Turn on **AutoFill Passwords and Passkeys**.
3. ☐ Under *AutoFill From*, turn **Locker** on. **TBD: should Passwords (Apple / iCloud Keychain) be turned OFF so Locker is the only option?**
✓ **Expected:** Locker shows a checkmark / is toggled on.

6.2 Locker — Live Activities

1. ☐ Go to **Settings** > **Apps** > **Locker** > **Live Activities** .
2. ☐ Turn **Live Activities** on. If *More Frequent Updates* appears, turn it on too. **TBD: confirm both should be ON**
✓ **Expected:** Live Activities toggle is green.

6.3 Cellular — 5G data mode

1. ☐ Go to **Settings** > **Cellular** > **Cellular Data Options** > **Data Mode** .
2. ☐ Select **Allow More Data on 5G**.
✓ **Expected:** Checkmark next to "Allow More Data on 5G".
This option only appears once cellular service is active. If it's missing, check the SIM/eSIM (Troubleshooting).

6.4 App Store — automatic downloads

1. ☐ Go to **Settings** > **Apps** > **App Store** .
2. ☐ Under *Automatic Downloads*: turn **App Downloads** on and **App Updates** on.

- ☐ Under *Cellular Data*: turn **Automatic Downloads** on, then tap **App Downloads** and choose **Always Allow**.

✓ **Expected:** Automatic Downloads on; Cellular → App Downloads = Always Allow.

7 Rover & Haiku — first launch

Each app only needs to be opened once so its permission prompts are answered. **Do not sign in.**

Rover

- ☐ Open **Rover**.
- ☐ Tap **Accept** / **Allow** on every prompt (notifications, camera, microphone, Bluetooth, local network, etc.).
- ☐ Stop at the sign-in screen. Swipe up to close the app.

Haiku

- ☐ Open **Haiku**.
- ☐ Tap **Accept** / **Allow** on every prompt.
- ☐ Stop at the sign-in screen. Swipe up to close the app.

TBD: if Rover or Haiku ask for a server / environment / organization code, what should CDW enter?

8 Quality check — every phone

Do this for each phone before it goes in a case. Record the result in the manifest. A phone that fails any check does not ship — set it aside and report it.

CHECK	WHERE	PASS WHEN
Label matches phone	Settings > General > About	Serial Number on screen = serial on the label
Managed by BCH	Settings > General > VPN & Device Management	A Boston Children's Hospital management profile is listed TBD: exact profile name
All apps installed	Home Screen	Every app on the list is present; none waiting or gray
Cellular service	Status bar; Settings > Cellular	Carrier name and signal bars show; phone number recorded TBD: test call or text required?
Hubb location	Settings > Apps > Hubb > Location	Always, Precise Location on
AutoFill	Settings > General > AutoFill & Passwords	Locker on
Live Activities	Settings > Apps > Locker	Live Activities on
5G data mode	Settings > Cellular > Cellular Data Options	Allow More Data on 5G
App Store downloads	Settings > Apps > App Store	Automatic Downloads on; Cellular App Downloads = Always Allow
Rover & Haiku	Open each app	Opens straight to the sign-in screen, no permission pop-ups
No Apple Account	Settings (top of screen)	Shows "Sign in to your iPhone" — nobody is signed in
Physical condition	Look at the phone	No scratches, cracks or dead pixels; buttons and cameras work

CHECK	WHERE	PASS WHEN
Battery	Status bar	At least TBD %

9 Case, pack & ship

Case

- ☐ Apply screen protector: **TBD: yes/no**
- ☐ Put the phone in its case. Make sure the label is flat and not wrinkled under the case.
- ☐ Power state for shipping: **TBD: power OFF, or leave ON**

Pack

- ☐ Retail boxes: **TBD: ship phones in their original Apple boxes, or discard/recycle boxes and bulk-pack?**
- ☐ Pack phones into **as few boxes as possible**, padded so they can't shift or rub screen-to-screen.
- ☐ Pack **all charging cables in a separate box** **TBD: power adapters too?**. Label it **CHARGING CABLES — QTY ____**.
- ☐ Number every box (**Box 1 of N**) and write the phone count on the outside.
- ☐ Put a **packing list** inside each box listing the serial numbers in that box. Enter the box number for each phone in the manifest.
- ☐ Email the completed manifest to BCH **before** the shipment leaves (see Contacts).
- ☐ Ship to the address BCH provided. Send the tracking number(s) to BCH.

A Device manifest

One row per phone, as an Excel or CSV file. BCH uses this to match every phone to its SIM, phone number and box.

COLUMN	EXAMPLE / NOTES
Serial Number	From Settings › General › About
IMEI	Settings › General › About
EID (eSIM) or ICCID (SIM card)	Whichever applies
Phone Number	Settings › Cellular › My Number
Carrier	
iOS Version	Settings › General › About
Enrolled in BCH management (Y/N)	From Quality check
All QA checks passed (Y/N)	If N, explain in Notes
Box Number	e.g. 3 of 12
Technician initials & date	
Notes	Damage, exceptions, anything unusual

B Troubleshooting & when to stop

STOP AND SET THE PHONE ASIDE WHEN...

- The Remote Management screen never appeared, or someone signed in with an Apple Account.
- The phone is damaged, won't power on, or has a serial not on BCH's order list.
- The same problem repeats after following the fix below once.

Record the serial number and the problem in the manifest Notes column and report it to BCH. Do not ship it.

PROBLEM	WHAT TO DO
No Remote Management screen — setup went straight to the Home Screen	The phone did not reach BCH's management. Erase it: Settings > General > Transfer or Reset iPhone > Erase All Content and Settings , then restart setup on a Wi-Fi network that meets the requirements. If it happens again, stop and report the serial — it may not be assigned in ABM.
"Unable to connect to the configuration server" or enrollment error	Check the Wi-Fi has internet and no sign-in page. Tap Try Again. If it fails twice, power off and report.
Apps stuck on Waiting... / Loading...	Confirm Wi-Fi is connected. Restart the phone (hold side + volume button, slide to power off, power on). Wait 15 more minutes. Still stuck: report the serial and app names.
An expected app never appears	Wait 30 minutes on Wi-Fi, restart once. Still missing: report — don't install it from the App Store.
No service / "No SIM" / SOS only	Physical SIM: power off, reseal the SIM. eSIM: check Settings > Cellular for a pending eSIM to activate. Still no service: record ICCID/EID and report.
"Allow More Data on 5G" is missing	Cellular service isn't active yet, or no 5G coverage at your site. Note it in the manifest; BCH will decide whether to ship.
Hubb location isn't "Always"	Settings > Apps > Hubb > Location → Always; turn on Precise Location.
Accidentally tapped "Don't Allow" in an app	Open Settings > Apps > (app name) and turn the permission back on.

PROBLEM**WHAT TO DO****Signed in with an Apple Account by mistake**

Stop. Erase the phone (see first row) and set it up again from Step 3. Tell BCH.

C Contacts**FOR****CONTACT**

Setup questions / problem phones

TBD: name, email, phone — BCH CMDI

Manifest & tracking numbers

TBD: email

Carrier / SIM activation issues

TBD: carrier account contact

Order / logistics (CDW side)

TBD: CDW project manager

Expected response time for questions: TBD. If BCH hasn't answered, set the affected phones aside and continue with the others.

D Per-phone bench checklist

Print one per phone (or per batch) and keep it with the device until it's packed.

Serial number _____ Phone number _____ ICCID / EID _____

Technician _____ Date _____ Box # _____

- | | |
|--|--|
| ☐ Serial checked against order list | ☐ Locker → Live Activities on |
| ☐ Charged $\geq$ 50% | ☐ Cellular → Allow More Data on 5G |
| ☐ Label printed & applied to phone back | ☐ App Store → Automatic Downloads on |
| ☐ SIM inserted / eSIM activated | ☐ App Store → Cellular App Downloads = Always Allow |
| ☐ English · United States | ☐ Rover opened — prompts allowed |
| ☐ Wi-Fi connected | ☐ Haiku opened — prompts allowed |
| ☐ Remote Management — enrolled | ☐ Label serial = Settings serial |
| ☐ No Apple Account signed in | ☐ BCH profile in VPN & Device Mgmt |
| ☐ Location Services enabled | ☐ Cellular service + number recorded |
| ☐ Hubb opened — permissions accepted | ☐ No damage |
| ☐ Hubb location = Always + Precise | ☐ Case on |
| ☐ All apps installed, none waiting | ☐ Manifest row complete |
| ☐ AutoFill → Locker on | ☐ Packed — box # recorded |